



Labcon, North America
3700 Lakeville Highway
Suite 200
Petaluma, CA 94954-5671 U.S.A.
707.766.2100
Fax 707.766.2199
www.labcon.com

Certificate of Analysis

Labcon North America certifies our low metal centrifuge tubes for use in applications requiring minimal trace metal conditions.

A 5% acetic acid solution was incubated in a test article for 18 hours at 55°C, followed by ICP analysis of trace elements found in the test solution.

The results indicate that the following elements fall below 1 ppb:

- *Calcium*
- *Magnesium*
- *Zinc*
- *Iron*
- *Manganese*
- *Copper*
- *Aluminum*
- *Silicon*
- *Nickel*
- *Vanadium*
- *Sodium*
- *Phosphorous*
- *Cobalt*
- *Chromium*
- *Potassium*
- *Lithium*
- *Lead*