



Now tested for Adenosine Triphosphate

Our 96 well, 0.2mL PCR* plates are available in styles to fit most popular thermocyclers. Autoclavable. Lot tested and certified free of Rnase, Dnase, DNA, and ATP. Polypropylene. See pages 88-89 for fit on your Thermocycler.



96-Well PCR* Plates
with printed Alpha Numeric Labeling

Item	Description	Qty/Pk	Sterile	Qty/Cs
3977-520	Natural	10		100



96-Well Low Profile PCR* Plates
with printed Alpha Numeric Labeling

Item	Description	Qty/Pk	Sterile	Qty/Cs
3971-520	Natural	20		100

(Low profile plates have smaller wells)

SEE COMPATIBILITY CHARTS
ON PAGES 88-89.




Barcoded plates are available incorporating
unique numbers.



96-Well Skirted PCR* Plates
with printed Alpha Numeric Labeling

Item	Description	Qty/Pk	Sterile	Qty/Cs
3978-520	Natural	10		100
3978-520-010	Natural with Barcode	10		100




96 High Well PCR* Plates
(these can be sealed with strip caps on page 81)

Item	Description	Qty/Pk	Sterile	Qty/Cs
3974-520	Natural	10		100

* PCR is a patented process of Hoffman-LaRoche. Use of the PCR process requires a license.

DNA tested


11



96-Well Half Skirt PCR* Plates
with printed Alpha Numeric Labeling

Item	Description	Qty/Pk	Sterile	Qty/Cs
3973-520	Natural	10		100
3973-520-010	Natural with Barcode	10		100

11



96-Well Low Profile Half Skirt PCR* Plates
with printed Alpha Numeric Labeling

Item	Description	Qty/Pk	Sterile	Qty/Cs
3976-520	Natural	10		100

(Low profile plates have smaller wells)

SEE COMPATIBILITY CHARTS
ON PAGES 88-89.


11



96-Well Half Skirt PCR* Plates, Raised Lip
with printed Alpha Numeric Labeling

Item	Description	Qty/Pk	Sterile	Qty/Cs
3979-520	Natural	10		100


19



96-Well Half Skirt PCR* Plates
for ABI style thermocyclers

Item	Description	Qty/Pk	Sterile	Qty/Cs
3972-520	Natural	10		100
3972-520-010	Natural with Barcode	10		100

20



96-Well PCR* Plates for Roche® Lightcycler® 480
with printed Alpha Numeric Labeling

Item	Description	Qty/Pk	Sterile	Qty/Cs
3962-520	Natural	10		100
3962-521	White	10		100

21



96-Well PCR* Plates for ABI® Fast 7500
with printed Alpha Numeric Labeling


Item	Description	Qty/Pk	Sterile	Qty/Cs
3961-520	Natural	10		100
3961-521	White	10		100

* PCR is a patented process of Hoffman-LaRoche. Use of the PCR process requires a license.

Now tested for Adenosine Triphosphate

384-Well PCR* Plates

384-well PCR* plates are made for high throughput thermocycling applications. These plates are flexible and skirted to allow bar coding. Each plate includes alpha numeric labeling with a frosted labeling area. Each well holds 25µL and has smooth inner walls for rapid heat transfer. Slightly raised walls seal with films or foils. Autoclavable. Made of polypropylene.



384-Well Skirted PCR* Plates
with Alpha Numeric Labeling

Item	Description	Qty/Pk	Sterile	Qty/Cs
3983-520	Natural	10		100
3983-520-010	Natural with Barcode	10		100

384-Well Skirted PCR* Plates for ABI Life Technologies Thermocyclers
with printed Alpha Numeric Labeling

Item	Description	Qty/Pk	Sterile	Qty/Cs
3984-520	Natural	10		100

SEE COMPATIBILITY
CHARTS
ON PAGES 88-89.


PCR Plate Sealing Mats

Sealing mats help reduce evaporation when cycling or storing for short periods. These reusable silicone mats fit most brands of PCR* plates and can be cleaned and reused. Autoclavable at 122°C.



96-Well Silicone Sealing Mats
Natural.

Item	Description	Qty/Pk	Sterile	Qty/Cs
3992-520	Round Wells	5		50



96-Well Silicone Sealing Mats
Natural.

Item	Description	Qty/Pk	Sterile	Qty/Cs
3994-520	Square Wells	5		50

* PCR is a patented process of Hoffman-LaRoche. Use of the PCR process requires a license.